

免疫质控品

货号: IA3111
批号: 2198EC

包装: 12 x 5 ml
效期: 2024-12

产品描述

本品适用于临床化学体外诊断的免疫质控，用于监测免疫测定的正确度或重复度。朗道供应三种浓度水平（水平 1: IA3109；水平 2: IA3110；水平 3: IA3111）。

安全预防措施和警告

本产品仅用于体外诊断。禁止用口吸。按照实验室常规预防措施对试剂进行处理。

该质控品采用人基质血清，对所有捐献者的血清均进行了 HIV(HIV1、HIV2)抗体、肝炎 B 表面抗原(HbsAg)和肝炎 C 病毒(HCV)抗体的测试，发现均呈阴性。所采用的方法均经 FDA 认证。

然而，既然没有一种方法能够完全保证其没有传染物质，因此该质控品和所有的病人样品均应当按照能够传播疾病的样品小心处理。

保存和稳定性

开瓶后，2~8℃保存。复溶后的血清如果保存在原始容器，并且无污染，可在 2~8℃保存 7 天，并可在-20℃保存 4 周，只能冻融一次。C 肽在 2~8℃可稳定 1 天。甲状旁腺素(PTH)和甲状腺球蛋白(TG)复溶后如果 2~8℃保存，要在 4 小时内检测；若-20℃保存，则须在 2 周内检测。促肾上腺皮质激素(ACTH)如果 2~8℃保存，要在 4 小时内检测。不提供促肾上腺皮质激素(ACTH)、醛固酮和 C-肽的冷冻稳定性说明。每次使用只吸取所需用量，剩余的样品不可返回原瓶。

未开瓶，2~8℃可保存至效期末。

复溶后的血清在受到细菌污染后，会导致许多组分稳定性的下降，应该将该瓶样品抛弃，重新复溶一瓶新的血清。

使用说明

该免疫质控品为冻干粉。

- 小心打开瓶盖，往瓶中准确加入 5 ml 蒸馏水，在 15~25℃的条件下复溶，盖上瓶盖，放置 30 分钟。轻轻翻转瓶子，确保内容物完全溶解；不要振荡，避免气泡的产生。
- 参照各个分析仪的质控说明。
- 冷藏未使用的材料，再次使用前，确保样品充分混合均匀。

需要自备的材料

移液管 蒸馏水

赋值

每一批质控血清都要送到参考实验室，这些实验室通过独特的统计分析方法得到一个结果进行赋值。质控范围值是平均值 ± 2 S.D.，质控结果应该落在这个区间内；然而，朗道赋值系统与客户检测系统存在的差异不可避免，所以靶值书仅供参考，建议每个实验室建立自己的参考靶值和可接受范围。

该质控品可用于醛固酮项目，但只提供水平 2 和水平 3 的靶值。

雅培 Abbott Architect 的 CA19-9 检验原理为抗原抗体反应，所使用的抗体材料是 1116-NA-19-9，有报道指出 1116-NA-19-9 的浓度很高时，Abbott Architect 系统的检验结果比其它方法的结果可能更高。

注：详细赋值信息请以原版英文说明书为准，原版说明书请在英国朗道公司官网 www.randox.com 进行下载。

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
17-羟基孕酮	nmol/l	39.2	29.4	49.0	Immunotech RIA
	ng/ml	13.0	9.73	16.3	
	nmol/l	49.2	36.9	61.5	ELISA
	ng/ml	16.3	12.2	20.4	
	nmol/l	42.2	31.7	52.8	Novatec ELISA
	ng/ml	14.0	10.5	17.5	
	nmol/l	53.6	40.2	67.0	DRG ELISA
	ng/ml	17.7	13.3	22.1	
	nmol/l	40.5	30.4	50.6	DIAsource RIA
	ng/ml	13.4	10.1	16.7	
	nmol/l	59.1	44.3	73.9	DIAsource ELISA
	ng/ml	19.6	14.7	24.5	
	nmol/l	57.9	43.4	72.4	LC/MS
	ng/ml	19.2	14.4	24.0	
25-OH-维生素D	nmol/l	27.9	20.9	34.9	IDS-Isys 17-OH Progesterone
	ng/ml	9.23	6.92	11.5	
	nmol/l	107	80.3	134	Diasorin Liaison
	nmol/l	140	105	175	Siemens Centaur XP/XPT/Classic
	nmol/l	94.2	70.7	118	Roche Vitamin D Total
	nmol/l	103	77.3	129	Abbott Architect (3L52)
	nmol/l	109	81.8	136	Vitros ECi
	nmol/l	84.8	63.6	106	Tosoh Series
	nmol/l	100	75.0	125	BioMerieux Vidas
	nmol/l	88.6	66.5	111	Beckman Dxl 600/800
	nmol/l	120	90.0	150	Fujirebio Lumipulse G Series
	nmol/l	90.3	67.7	113	Beckman Access 25 OH Vitamin D Total
	nmol/l	111	83.3	139	LC/MS
	nmol/l	102	76.5	128	Abbott Architect (5P02)/Alinity (8P45)
	nmol/l	107	80.3	134	Diasorin Liaison XL
	nmol/l	107	80.3	134	Roche Vitamin D Total II
	nmol/l	107	80.3	134	Abbott Alinity i
	nmol/l	143	107	179	Siemens Atellica IM
	nmol/l	102	76.5	128	Roche Vitamin D Total II e801
	促肾上腺皮质激素 (ACTH)	pmol/l	49.0	36.8	61.3
pg/ml		223	167	279	
pmol/l		28.1	21.1	35.1	Roche Cobas 4000/E411
pg/ml		128	95.9	160	
pmol/l		52.0	39.0	65.0	Siemens Immulite 1000
pg/ml		236	177	295	
pmol/l		29.9	22.4	37.4	Roche Cobas e402/e801
pg/ml		136	102	170	
pmol/l		16.7	12.5	20.9	Tosoh AIA-CL Series
pg/ml		75.9	56.8	95.0	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学	
醛固酮	pmol/l	1111	833	1389	DRG ELISA	
	pg/ml	401	301	501		
	pmol/l	265	199	331	CIS RIA	
	pg/ml	95.6	71.8	119		
	pmol/l	814	611	1018	IDS iSYS	
	pg/ml	294	220	368		
	pmol/l	667	500	834	Diasorin Liaison	
	pg/ml	241	180	302		
	pmol/l	688	516	860	Diasorin Liaison XL	
	pg/ml	248	186	310		
	甲胎蛋白	KIU/l = IU/ml	167	134	200	Abbott AFP-EIA
		ng/ml	202	162	242	
KIU/l = IU/ml		166	133	199	Abbott Architect	
ng/ml		201	161	241		
KIU/l = IU/ml		188	150	226	BioMerieux Vidas	
ng/ml		227	182	272		
KIU/l = IU/ml		192	154	230	Siemens Centaur XP/XPT/Classic	
ng/ml		232	186	278		
KIU/l = IU/ml		172	138	206	Siemens Immulite 2000/2500	
ng/ml		208	167	249		
KIU/l = IU/ml		184	147	221	Siemens Immulite 1000	
ng/ml		223	178	268		
KIU/l = IU/ml		158	126	190	Beckman DxI800	
ng/ml		191	152	230		
KIU/l = IU/ml		181	145	217	Roche Elecsys	
ng/ml		219	175	263		
KIU/l = IU/ml		173	138	208	Diasorin Liaison	
ng/ml		209	167	251		
KIU/l = IU/ml		165	132	198	Beckman Access	
ng/ml		200	160	240		
KIU/l = IU/ml		157	126	188	Tosoh Series	
ng/ml		190	152	228		
KIU/l = IU/ml		176	141	211	Vitros ECi	
ng/ml		183	146	220		
KIU/l = IU/ml		187	150	224	Roche Cobas e601/602	
ng/ml		226	182	270		
KIU/l = IU/ml		263	210	316	CIS IRMA	
ng/ml		318	254	382		
KIU/l = IU/ml		189	151	227	Roche Cobas 4000/E411	
ng/ml		229	183	275		
KIU/l = IU/ml	181	145	217	Fujirebio Lumipulse G Series		
ng/ml	219	175	263			
KIU/l = IU/ml	177	142	212	Siemens Dimension		
ng/ml	214	172	256			
KIU/l = IU/ml	185	148	222	Siemens Centaur CP		
ng/ml	224	179	269			

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学	
甲胎蛋白	KIU/l = IU/ml	199	159	239	Sysmex HISCL Series	
	ng/ml	241	192	290		
	KIU/l = IU/ml	183	146	220	Diasorin Liaison XL	
	ng/ml	221	177	265		
	KIU/l = IU/ml	188	150	226	Siemens Atellica IM	
	ng/ml	227	182	272		
	KIU/l = IU/ml	165	132	198	Abbott Alinity i	
	ng/ml	200	160	240		
	KIU/l = IU/ml	179	143	215	Roche Cobas e402/e801	
	ng/ml	217	173	261		
	KIU/l = IU/ml	158	126	190	Tosoh AIA-CL Series	
	ng/ml	191	152	230		
丁胺卡那霉素	μmol/l	21.8	17.4	26.2	EMIT	
	μg/ml	12.8	10.2	15.4		
	μmol/l	23.1	18.5	27.7	Roche Cobas 6000/8000	
	μg/ml	13.5	10.8	16.2		
	μmol/l	24.5	19.6	29.4	Abbott Architect	
	μg/ml	14.3	11.5	17.1		
	μmol/l	24.0	19.2	28.8	Thermo Scientific QMS Amikacin	
	μg/ml	14.1	11.2	17.0		
	μmol/l	24.9	19.9	29.9	Abbott Alinity	
	μg/ml	14.6	11.7	17.5		
	μmol/l	24.1	19.3	28.9	Gravimetric	
	μg/ml	14.1	11.3	16.9		
雄烯二酮	nmol/l	24.0	18.0	30.0	Beckman Coulter RIA	
	ng/ml	6.89	5.17	8.61		
	nmol/l	23.3	17.5	29.1	DRG ELISA	
	ng/ml	6.69	5.02	8.36		
	nmol/l	24.1	18.1	30.1	LC/MS	
	ng/ml	6.92	5.19	8.65		
	nmol/l	32.3	24.2	40.4	Diasorin Liaison XL	
	ng/ml	9.27	6.95	11.6		
	β-2-微球蛋白	μg/ml = mg/l	6.81	5.45	8.17	Abbott Architect
		μg/ml = mg/l	6.93	5.54	8.32	BioMerieux Vidas
μg/ml = mg/l		6.72	5.38	8.06	Behring Nephelometer	
μg/ml = mg/l		5.76	4.61	6.91	Siemens Immulite 2000/2500	
μg/ml = mg/l		6.61	5.29	7.93	Tosoh Series	
μg/ml = mg/l		6.46	5.17	7.75	Roche Cobas 6000/8000	
μg/ml = mg/l		5.55	4.44	6.66	Fujirebio Lumipulse G Series	
μg/ml = mg/l		6.23	4.98	7.48	Diasorin Liaison XL	
μg/ml = mg/l		7.43	5.94	8.92	Siemens Atellica CH	
μg/ml = mg/l		6.80	5.44	8.16	Abbott Alinity	
C-肽	nmol/l	3.390	2.540	4.240	Siemens Immulite 2000/2500	
	ng/ml	10.2	7.67	12.7		
	nmol/l	3.220	2.420	4.030	Roche Elecsys	
	ng/ml	9.72	7.31	12.1		

免疫质控品-水平3

货号: IA3111 批号: 2198EC 规格: 12 x 5ml 效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
C-肽	nmol/l	3.260	2.450	4.080	Roche Cobas e601/602
	ng/ml	9.84	7.40	12.3	
	nmol/l	3.710	2.780	4.640	Siemens Immulite 1000
	ng/ml	11.2	8.39	14.0	
	nmol/l	3.220	2.420	4.030	Diasorin Liaison
	ng/ml	9.72	7.31	12.1	
	nmol/l	3.290	2.470	4.110	Roche Cobas 4000/E411
	ng/ml	9.93	7.46	12.4	
	nmol/l	3.070	2.300	3.840	Abbott Architect
	ng/ml	9.27	6.94	11.6	
	nmol/l	2.370	1.780	2.960	Siemens Centaur XP/XPT/Classic
	ng/ml	7.16	5.37	8.95	
	nmol/l	2.860	2.150	3.580	Tosoh Series
	ng/ml	8.63	6.49	10.8	
	nmol/l	3.460	2.600	4.330	Diasorin Liaison XL
	ng/ml	10.4	7.85	13.0	
	nmol/l	3.180	2.390	3.980	Roche Cobas e402/e801
	ng/ml	9.60	7.22	12.0	
	nmol/l	2.290	1.720	2.860	Siemens Atellica IM
	ng/ml	6.91	5.19	8.63	
nmol/l	3.050	2.290	3.810	Vitros ECi	
ng/ml	9.21	6.91	11.5		
nmol/l	3.530	2.650	4.410	Tosoh AIA-CL Series	
ng/ml	10.7	8.00	13.4		
nmol/l	2.120	1.590	2.650	Beckman Access	
ng/ml	6.40	4.80	8.00		
nmol/l	2.900	2.180	3.630	Abbott Alinity i	
ng/ml	8.76	6.58	10.9		
CA 15-3	U/ml	114	91.2	137	Abbott Architect
	U/ml	68.1	54.5	81.7	BioMerieux Vidas
	U/ml	76.0	60.8	91.2	Siemens Centaur XP/XPT/Classic
	U/ml	83.9	67.1	101	Siemens Immulite 2000/2500
	U/ml	103	82.4	124	Siemens Immulite 1000
	U/ml	33.1	26.5	39.7	Beckman Dxl800
	U/ml	75.9	60.7	91.1	Roche Elecsys
	U/ml	82.5	66.0	99.0	Diasorin Liaison
	U/ml	32.4	25.9	38.9	Beckman Access
	U/ml	94.2	75.4	113	Tosoh Series
	U/ml	67.9	54.3	81.5	Vitros ECi
	U/ml	74.9	59.9	89.9	Roche Cobas e601/602
	U/ml	75.0	60.0	90.0	Roche Cobas 4000/E411
	U/ml	88.3	70.6	106	Fujirebio Lumipulse G Series
	U/ml	78.6	62.9	94.3	Siemens Centaur CP
	U/ml	87.7	70.2	105	Diasorin Liaison XL
	U/ml	84.2	67.4	101	Siemens Atellica IM

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
CA 15-3	U/ml	71.6	57.3	85.9	Roche Cobas e402/e801
	U/ml	108	86.4	130	Abbott Alinity i
CA 19-9	U/ml	285	228	342	BioMerieux Vidas
	U/ml	339	271	407	Siemens Centaur XP/XPT/Classic
	U/ml	361	289	433	Siemens Immulite 2000/2500
	U/ml	347	278	416	Beckman Dxl800
	U/ml	160	128	192	Roche Elecsys
	U/ml	311	249	373	Diasorin Liaison
	U/ml	388	310	466	Beckman Access
	U/ml	130	104	156	Tosoh Series
	U/ml	354	283	425	Vitros ECI
	U/ml	519	415	623	Abbott Architect XR/XR2/XRL
	U/ml	162	130	194	Roche Cobas e601/602
	U/ml	164	131	197	Roche Cobas 4000/E411
	U/ml	347	278	416	Fujirebio Lumipulse G Series
	U/ml	206	165	247	Izotop RIA
	U/ml	261	209	313	Siemens Centaur CP
	U/ml	302	242	362	Diasorin Liaison XL
	U/ml	340	272	408	Siemens Atellica IM
	U/ml	162	130	194	Roche Cobas e402/e801
	U/ml	528	422	634	Abbott Alinity XR/XR2/XRL
CA125	U/ml	241	193	289	Abbott Architect
	U/ml	197	158	236	BioMerieux Vidas
	U/ml	201	161	241	Siemens Centaur XP/XPT/Classic OV Kit
	U/ml	206	165	247	Siemens Centaur XP/XPT/Classic CA125 II kit
	U/ml	176	141	211	Siemens Immulite 2000/2500
	U/ml	208	166	250	Siemens Immulite 1000
	U/ml	158	126	190	Beckman Dxl800
	U/ml	176	141	211	Roche Elecsys
	U/ml	157	126	188	Beckman Access
	U/ml	330	264	396	Tosoh Series
	U/ml	188	150	226	Vitros ECI
	U/ml	182	146	218	Roche Cobas 4000/E411
	U/ml	183	146	220	Roche Cobas e601/602
	U/ml	159	127	191	Fujirebio Lumipulse G Series
	U/ml	199	159	239	Siemens Centaur CP CA125 II kit
	U/ml	199	159	239	Siemens Centaur CP OV kit
	U/ml	242	194	290	Diasorin Liaison XL
	U/ml	232	186	278	Siemens Atellica IM
	U/ml	175	140	210	Roche Cobas e402/e801
	U/ml	247	198	296	Abbott Alinity i
立痛定 (卡马西平)	μmol/l	71.0	56.8	85.2	Abbott Architect (CMIA)
	μg/ml	16.8	13.4	20.2	
	μmol/l	72.8	58.2	87.4	Siemens Centaur XP/XPT/Classic
	μg/ml	17.2	13.8	20.6	
	μmol/l	64.9	51.9	77.9	Siemens Dimension
μg/ml	15.3	12.3	18.3		

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
立痛定 (卡马西平)	μmol/l	59.9	47.9	71.9	Siemens Immulite 2000/2500
	μg/ml	14.2	11.3	17.1	
	μmol/l	73.7	59.0	88.4	Siemens Immulite 1000
	μg/ml	17.4	13.9	20.9	
	μmol/l	74.6	59.7	89.5	Roche Integra
	μg/ml	17.6	14.1	21.1	
	μmol/l	67.3	53.8	80.8	Vitros
	μg/ml	15.9	12.7	19.1	
	μmol/l	70.6	56.5	84.7	Roche Cobas 6000/8000
	μg/ml	16.7	13.4	20.0	
	μmol/l	71.1	56.9	85.3	Beckman DxC
	μg/ml	16.8	13.5	20.1	
	μmol/l	74.2	59.4	89.0	Beckman AU400/600/640/2700/5400
	μg/ml	17.5	14.0	21.0	
	μmol/l	56.3	45.0	67.6	Abbott Architect (PETINA)
	μg/ml	13.3	10.6	16.0	
	μmol/l	79.5	63.6	95.4	EMIT
	μg/ml	18.8	15.0	22.6	
	μmol/l	67.0	53.6	80.4	Siemens Atellica CH
	μg/ml	15.8	12.7	18.9	
μmol/l	49.5	39.6	59.4	Abbott Alinity (PETINA)	
μg/ml	11.7	9.36	14.0		
μmol/l	72.2	57.8	86.6	Gravimetric	
μg/ml	17.1	13.7	20.5		
癌胚抗原(CEA)	ng/ml = μg/l	77.6	62.1	93.1	Abbott Architect
	ng/ml = μg/l	58.3	46.6	70.0	BioMerieux Vidas
	ng/ml = μg/l	54.5	43.6	65.4	Siemens Centaur XP/XPT/Classic
	ng/ml = μg/l	55.0	44.0	66.0	CIS RIA Coated tube
	ng/ml = μg/l	81.1	64.9	97.3	Siemens Immulite 2000/2500
	ng/ml = μg/l	81.5	65.2	97.8	Siemens Immulite 1000
	ng/ml = μg/l	53.8	43.0	64.6	Beckman Dxl800
	ng/ml = μg/l	52.4	41.9	62.9	Roche Elecsys
	ng/ml = μg/l	52.2	41.8	62.6	Diasorin Liaison
	ng/ml = μg/l	57.0	45.6	68.4	Beckman Access
	ng/ml = μg/l	86.6	69.3	104	Tosoh Series
	ng/ml = μg/l	65.8	52.6	79.0	Vitros ECI
	ng/ml = μg/l	51.6	41.3	61.9	Roche Cobas e601/602
	ng/ml = μg/l	51.9	41.5	62.3	Roche Cobas 4000/E411
	ng/ml = μg/l	72.7	58.2	87.2	Fujirebio Lumipulse G Series
	ng/ml = μg/l	54.4	43.5	65.3	Siemens Centaur CP
	ng/ml = μg/l	55.4	44.3	66.5	Diasorin Liaison XL
	ng/ml = μg/l	53.4	42.7	64.1	Siemens Atellica IM
	ng/ml = μg/l	52.5	42.0	63.0	Roche Cobas e402/e801
	ng/ml = μg/l	76.5	61.2	91.8	Abbott Alinity i
皮质醇	nmol/l	861	646	1076	Abbott Architect
	μg/dl	31.0	23.3	38.7	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
皮质醇	nmol/l	1010	758	1263	BioMerieux Vidas
	µg/dl	36.4	27.3	45.5	
	nmol/l	1093	820	1366	Siemens Centaur XP/XPT/Classic
	µg/dl	39.3	29.5	49.1	
	nmol/l	789	592	986	CIS RIA
	µg/dl	28.4	21.3	35.5	
	nmol/l	863	647	1079	Siemens Immulite 2000/2500
	µg/dl	31.1	23.3	38.9	
	nmol/l	908	681	1135	Siemens Immulite 1000
	µg/dl	32.7	24.5	40.9	
	nmol/l	984	738	1230	Beckman Dxl800
	µg/dl	35.4	26.6	44.2	
	nmol/l	900	675	1125	Roche Elecsys
	µg/dl	32.4	24.3	40.5	
	nmol/l	1058	794	1323	Beckman Access
	µg/dl	38.1	28.6	47.6	
	nmol/l	1042	782	1303	Tosoh Series
	µg/dl	37.5	28.2	46.8	
	nmol/l	898	674	1123	Vitros ECi
	µg/dl	32.3	24.3	40.3	
nmol/l	961	721	1201	Roche Cobas e601/602	
µg/dl	34.6	26.0	43.2		
nmol/l	990	743	1238	Roche Cobas 4000/E411	
µg/dl	35.6	26.7	44.5		
nmol/l	1192	894	1490	Siemens Centaur CP	
µg/dl	42.9	32.2	53.6		
nmol/l	988	741	1235	Roche Cobas 4000/e411 Gen II	
µg/dl	35.6	26.7	44.5		
nmol/l	945	709	1181	Roche Cobas e601/602 Gen II	
µg/dl	34.0	25.5	42.5		
nmol/l	955	716	1194	Roche Elecsys Gen II	
µg/dl	34.4	25.8	43.0		
nmol/l	1038	779	1298	Diasorin Liaison XL	
µg/dl	37.4	28.0	46.8		
nmol/l	1158	869	1448	Siemens Atellica IM	
µg/dl	41.7	31.3	52.1		
nmol/l	966	725	1208	Roche Cobas e402/e801	
µg/dl	34.8	26.1	43.5		
nmol/l	1092	819	1365	Tosoh AIA-CL Series	
µg/dl	39.3	29.5	49.1		
nmol/l	847	635	1059	Abbott Alinity i	
µg/dl	30.5	22.9	38.1		
去表雄酮硫酸盐 (DHEA-S)	µmol/l	19.4	15.5	23.3	Abbott Architect
	µmol/l	16.7	13.4	20.0	Siemens Immulite 2000/2500
	µmol/l	15.7	12.6	18.8	Siemens Immulite 1000

免疫质控品-水平3

货号: IA3111 批号: 2198EC 规格: 12 x 5ml 效期: 2024-12-28

范围

分析物	单位	靶值	低值	高值	方法学
去表雄酮硫酸盐 (DHEA-S)	μmol/l	15.4	12.3	18.5	Beckman DxI800
	μmol/l	25.7	20.6	30.8	Roche Elecsys
	μmol/l	16.8	13.4	20.2	Beckman Access
	μmol/l	15.3	12.2	18.4	Diasorin Liaison
	μmol/l	14.7	11.8	17.6	Siemens Centaur XP/XPT/Classic
	μmol/l	24.8	19.8	29.8	Roche Cobas 4000/E411
	μmol/l	24.3	19.4	29.2	Roche Cobas e601/602
	μmol/l	16.6	13.3	19.9	LC/MS
	μmol/l	14.8	11.8	17.8	Siemens Atellica IM
	μmol/l	23.3	18.6	28.0	Roche Cobas e402/e801
μmol/l	19.8	15.8	23.8	Abbott Alinity i	
地高辛	nmol/l	3.98	3.18	4.78	Abbott Architect (CMIA)
	ng/ml	3.11	2.48	3.74	
	nmol/l	4.39	3.51	5.27	BioMerieux Vidas
	ng/ml	3.43	2.74	4.12	
	nmol/l	5.35	4.28	6.42	Siemens Centaur XP/XPT/Classic
	ng/ml	4.18	3.34	5.02	
	nmol/l	4.66	3.73	5.59	Siemens Dimension
	ng/ml	3.64	2.91	4.37	
	nmol/l	5.21	4.17	6.25	Siemens Immulite 2000/2500
	ng/ml	4.07	3.26	4.88	
	nmol/l	4.70	3.76	5.64	Beckman DxI800
	ng/ml	3.67	2.94	4.40	
	nmol/l	4.43	3.54	5.32	Roche Integra
	ng/ml	3.46	2.76	4.16	
	nmol/l	4.68	3.74	5.62	Roche Cobas e601/602
	ng/ml	3.66	2.92	4.40	
	nmol/l	5.09	4.07	6.11	Beckman AU400/600/640/2700/5400
	ng/ml	3.98	3.18	4.78	
	nmol/l	4.90	3.92	5.88	Roche Cobas 4000/E411
	ng/ml	3.83	3.06	4.60	
nmol/l	4.38	3.50	5.26	Abbott Architect (PETINA)	
ng/ml	3.42	2.73	4.11		
nmol/l	4.67	3.74	5.60	Roche Cobas e402/e801	
ng/ml	3.65	2.92	4.38		
nmol/l	5.78	4.62	6.94	Siemens Atellica IM	
ng/ml	4.51	3.61	5.41		
nmol/l	4.53	3.62	5.44	Siemens Atellica CH	
ng/ml	3.54	2.83	4.25		
nmol/l	4.55	3.64	5.46	Abbott Alinity (CMIA)	
ng/ml	3.55	2.84	4.26		
nmol/l	4.55	3.64	5.46	Abbott Alinity (PETINA)	
ng/ml	3.55	2.84	4.26		
nmol/l	5.12	4.10	6.14	Gravimetric	
ng/ml	4.00	3.20	4.80		
乙琥胺	μmol/l	493	394	592	Gravimetric
	μg/ml	69.9	55.9	83.9	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
血清铁蛋白	ng/ml = µg/l	338	270	406	Abbott Architect
	ng/ml = µg/l	244	195	293	BioMerieux Vidas
	ng/ml = µg/l	239	191	287	Siemens Centaur XP/XPT/Classic
	ng/ml = µg/l	289	231	347	Siemens Dimension
	ng/ml = µg/l	266	213	319	Siemens Immulite 2000/2500
	ng/ml = µg/l	228	182	274	Siemens Immulite 1000
	ng/ml = µg/l	214	171	257	Beckman Dxl800
	ng/ml = µg/l	315	252	378	Roche Elecsys/Cobas e411
	ng/ml = µg/l	301	241	361	Diasorin Liaison
	ng/ml = µg/l	310	248	372	Roche Modular E170 / Roche Cobas e601/602
	ng/ml = µg/l	221	177	265	Beckman Access
	ng/ml = µg/l	221	177	265	Tosoh Series
	ng/ml = µg/l	192	154	230	Vitros ECi
	ng/ml = µg/l	307	246	368	Roche Cobas 4000/c311
	ng/ml = µg/l	250	200	300	Beckman AU400/600/640/2700/5400
	ng/ml = µg/l	295	236	354	Roche Cobas c501/502/503/701/702
	ng/ml = µg/l	233	186	280	Siemens Centaur CP
	ng/ml = µg/l	326	261	391	Diasorin Liaison XL
	ng/ml = µg/l	346	277	415	Abbott Alinity i
ng/ml = µg/l	230	184	276	Siemens Atellica IM	
ng/ml = µg/l	309	247	371	Roche Cobas e402/e801	
叶酸	nmol/l	17.9	13.6	22.2	Abbott Architect
	ng/ml	7.89	6.00	9.78	
	nmol/l	39.6	30.1	49.1	Siemens Centaur XP/XPT/Classic (enhanced kit)
	ng/ml	17.5	13.3	21.7	
	nmol/l	15.0	11.4	18.6	Siemens Immulite 2000
	ng/ml	6.62	5.03	8.21	
	nmol/l	21.6	16.4	26.8	Beckman Dxl 600 / 800 - A14208
	ng/ml	9.53	7.23	11.8	
	nmol/l	21.7	16.5	26.9	Beckman Access
	ng/ml	9.57	7.28	11.9	
	nmol/l	29.5	22.4	36.6	Tosoh Series
	ng/ml	13.0	9.88	16.1	
	nmol/l	14.0	10.6	17.4	Vitros ECi
	ng/ml	6.17	4.67	7.67	
	nmol/l	12.2	9.27	15.1	Siemens Dimension
	ng/ml	5.38	4.09	6.67	
	nmol/l	21.6	16.4	26.8	Beckman Dxl 600 / 800 - A98032
	ng/ml	9.53	7.23	11.8	
	nmol/l	44.0	33.4	54.6	Siemens Centaur CP (enhanced kit)
	ng/ml	19.4	14.7	24.1	
nmol/l	22.4	17.0	27.8	Roche Folate III - 04476433	
ng/ml	9.88	7.50	12.3		
nmol/l	22.9	17.4	28.4	Roche Elecsys Folate	
ng/ml	10.1	7.67	12.5		

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学	
叶酸	nmol/l	22.9	17.4	28.4	Roche Folate III - 07559992	
	ng/ml	10.1	7.67	12.5		
	nmol/l	34.8	26.4	43.2	Siemens Atellica IM	
	ng/ml	15.3	11.6	19.0		
	nmol/l	19.0	14.4	23.6	Abbott Alinity i	
	ng/ml	8.38	6.35	10.4		
	nmol/l	23.0	17.5	28.5	Roche Folate 07027290 e801	
	ng/ml	10.1	7.72	12.5		
	游离 T3	pmol/l	20.7	15.5	25.9	Abbott Architect, 2 point calibration
		ng/dl	1.35	1.01	1.69	
pg/ml		13.5	10.1	16.9		
pmol/l		21.1	15.8	26.4	BioMerieux Vidas	
ng/dl		1.37	1.03	1.71		
pg/ml		13.7	10.3	17.1		
pmol/l		20.8	15.6	26.0	Siemens Centaur XP/XPT/Classic	
ng/dl		1.35	1.01	1.69		
pg/ml		13.5	10.2	16.8		
pmol/l		13.4	10.1	16.8	Siemens Immulite 2000/2500	
ng/dl		0.871	0.657	1.09		
pg/ml		8.72	6.58	10.9		
pmol/l		13.4	10.1	16.8	Siemens Immulite 1000	
ng/dl		0.871	0.657	1.09		
pg/ml		8.72	6.58	10.9		
pmol/l		15.1	11.3	18.9	Beckman Dxl800	
ng/dl		0.982	0.735	1.23		
pg/ml		9.83	7.36	12.3		
pmol/l		26.3	19.7	32.9	Roche Elecsys	
ng/dl		1.71	1.28	2.14		
pg/ml		17.1	12.8	21.4		
pmol/l		28.0	21.0	35.0	Diasorin Liaison	
ng/dl		1.82	1.37	2.27		
pg/ml		18.2	13.7	22.7		
pmol/l		16.6	12.5	20.8	Beckman Access	
ng/dl		1.08	0.813	1.35		
pg/ml		10.8	8.14	13.5		
pmol/l		33.5	25.1	41.9	Tosoh Series	
ng/dl		2.18	1.63	2.73		
pg/ml		21.8	16.3	27.3		
pmol/l	35.2	26.4	44.0	Vitros ECi		
ng/dl	2.29	1.72	2.86			
pg/ml	22.9	17.2	28.6			
pmol/l	26.0	19.5	32.5	Roche Cobas 4000/E411		
ng/dl	1.69	1.27	2.11			
pg/ml	16.9	12.7	21.1			
pmol/l	26.2	19.7	32.8	Roche Cobas e601/602		
ng/dl	1.70	1.28	2.12			
pg/ml	17.1	12.8	21.4			

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
游离 T3	pmol/l	26.1	19.6	32.6	Fujirebio Lumipulse G Series
	ng/dl	1.70	1.27	2.13	
	pg/ml	17.0	12.8	21.2	
	pmol/l	20.7	15.5	25.9	Siemens Centaur CP
	ng/dl	1.35	1.01	1.69	
	pg/ml	13.5	10.1	16.9	
	pmol/l	24.4	18.3	30.5	Sysmex HISCL Series
	ng/dl	1.59	1.19	1.99	
	pg/ml	15.9	11.9	19.9	
	pmol/l	31.6	23.7	39.5	Diasorin Liaison XL
	ng/dl	2.05	1.54	2.56	
	pg/ml	20.6	15.4	25.8	
	pmol/l	22.6	17.0	28.3	Siemens Atellica IM
	ng/dl	1.47	1.11	1.83	
	pg/ml	14.7	11.1	18.3	
	pmol/l	21.5	16.1	26.9	Abbott Alinity, 2 point cal
	ng/dl	1.40	1.05	1.75	
	pg/ml	14.0	10.5	17.5	
	pmol/l	21.1	15.8	26.4	Abbott Architect, 6 point calibration
	ng/dl	1.37	1.03	1.71	
	pg/ml	13.7	10.3	17.1	
	pmol/l	26.1	19.6	32.6	Roche Cobas e402/e801
	ng/dl	1.70	1.27	2.13	
	pg/ml	17.0	12.8	21.2	
pmol/l	25.4	19.1	31.8	Tosoh AIA-CL Series	
ng/dl	1.65	1.24	2.06		
pg/ml	16.5	12.4	20.6		
pmol/l	21.6	16.2	27.0	Abbott Alinity, 6 point cal	
ng/dl	1.40	1.05	1.75		
pg/ml	14.1	10.5	17.7		
游离 T4	pmol/l	45.0	33.8	56.3	Abbott Architect
	ng/dl	3.51	2.64	4.38	
	pg/ml	35.1	26.4	43.8	
	pmol/l	46.2	34.7	57.8	Siemens Centaur XP/XPT/Classic
	ng/dl	3.60	2.71	4.49	
	pg/ml	36.0	27.1	44.9	
	pmol/l	54.4	40.8	68.0	CIS RIA Coated tube
	ng/dl	4.24	3.18	5.30	
	pg/ml	42.4	31.8	53.0	
	pmol/l	69.5	52.1	86.9	Siemens Immulite 2000/2500
	ng/dl	5.42	4.06	6.78	
	pg/ml	54.2	40.6	67.8	
	pmol/l	72.5	54.4	90.6	Siemens Immulite 1000
	ng/dl	5.66	4.24	7.08	
	pg/ml	56.6	42.4	70.8	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
游离 T4	pmol/l	53.9	40.4	67.4	Beckman Dxl800
	ng/dl	4.20	3.15	5.25	
	pg/ml	42.0	31.5	52.5	
	pmol/l	63.0	47.3	78.8	Roche Elecsys
	ng/dl	4.91	3.69	6.13	
	pg/ml	49.1	36.9	61.3	
	pmol/l	80.9	60.7	101	Diasorin Liaison
	ng/dl	6.31	4.73	7.89	
	pg/ml	63.1	47.3	78.9	
	pmol/l	49.7	37.3	62.1	Beckman Access
	ng/dl	3.88	2.91	4.85	
	pg/ml	38.8	29.1	48.5	
	pmol/l	67.2	50.4	84.0	Tosoh Series
	ng/dl	5.24	3.93	6.55	
	pg/ml	52.4	39.3	65.5	
	pmol/l	89.8	67.4	112	Vitros ECi
	ng/dl	7.00	5.26	8.74	
	pg/ml	70.0	52.6	87.4	
	pmol/l	65.6	49.2	82.0	Roche Cobas e601/602
	ng/dl	5.12	3.84	6.40	
	pg/ml	51.2	38.4	64.0	
	pmol/l	48.3	36.2	60.4	Immunotech IRMA
	ng/dl	3.77	2.82	4.72	
	pg/ml	37.7	28.2	47.2	
	pmol/l	63.2	47.4	79.0	Roche Cobas 4000/E411
	ng/dl	4.93	3.70	6.16	
	pg/ml	49.3	37.0	61.6	
	pmol/l	65.9	49.4	82.4	Biomerieux Vidas FT4N Kit
	ng/dl	5.14	3.85	6.43	
	pg/ml	51.4	38.5	64.3	
	pmol/l	61.3	46.0	76.6	Fujirebio Lumipulse G Series
	ng/dl	4.78	3.59	5.97	
	pg/ml	47.8	35.9	59.7	
	pmol/l	49.2	36.9	61.5	Siemens Centaur CP
	ng/dl	3.84	2.88	4.80	
	pg/ml	38.4	28.8	48.0	
	pmol/l	75.4	56.6	94.3	Sysmex HISCL Series
	ng/dl	5.88	4.41	7.35	
	pg/ml	58.8	44.1	73.5	
	pmol/l	99.0	74.3	124	Diasorin Liaison XL
	ng/dl	7.72	5.80	9.64	
	pg/ml	77.2	58.0	96.4	
	pmol/l	49.4	37.1	61.8	Siemens Atellica IM
	ng/dl	3.85	2.89	4.81	
	pg/ml	38.5	28.9	48.1	
	pmol/l	41.4	31.1	51.8	Abbott Alinity i
	ng/dl	3.23	2.43	4.03	
	pg/ml	32.3	24.3	40.3	

RANDOX

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
游离T4	pmol/l	65.7	49.3	82.1	Roche Cobas e402/e801
	ng/dl	5.12	3.85	6.39	
	pg/ml	51.2	38.5	63.9	
	pmol/l	60.1	45.1	75.1	Tosoh AIA-CL Series
	ng/dl	4.69	3.52	5.86	
	pg/ml	46.9	35.2	58.6	
卵泡刺激素 (FSH)	mU/ml	54.1	43.3	64.9	Abbott Architect
	mU/ml	70.5	56.4	84.6	BioMerieux Vidas
	mU/ml	81.0	64.8	97.2	Siemens Centaur XP/XPT/Classic
	mU/ml	62.6	50.1	75.1	Siemens Immulite 2000/2500
	mU/ml	64.4	51.5	77.3	Siemens Immulite 1000
	mU/ml	65.5	52.4	78.6	Beckman Dxl800
	mU/ml	59.2	47.4	71.0	Roche/BM Elecsys
	mU/ml	77.9	62.3	93.5	Diasorin Liaison
	mU/ml	64.2	51.4	77.0	Beckman Access
	mU/ml	59.6	47.7	71.5	Tosoh Series
	mU/ml	45.5	36.4	54.6	Vitros Eci
	mU/ml	61.0	48.8	73.2	Roche Cobas 4000/E411
	mU/ml	55.3	44.2	66.4	Roche Cobas e601/602
	mU/ml	52.7	42.2	63.2	Fujirebio Lumipulse G Series
	mU/ml	71.6	57.3	85.9	Siemens Centaur CP
	mU/ml	89.6	71.7	108	Diasorin Liaison XL
	mU/ml	80.3	64.2	96.4	Siemens Atellica IM
	mU/ml	54.2	43.4	65.0	Abbott Alinity i
	mU/ml	56.3	45.0	67.6	Roche Cobas e402/e801
庆大霉素	μmol/l	10.6	8.48	12.7	Roche Integra
	μg/ml	5.07	4.05	6.09	
	μmol/l	13.7	11.0	16.4	Siemens Centaur XP/XPT/Classic
	μg/ml	6.55	5.26	7.84	
	μmol/l	11.5	9.20	13.8	Abbott Architect (CMIA)
	μg/ml	5.50	4.40	6.60	
	μmol/l	12.4	9.92	14.9	Siemens Dimension Vista
	μg/ml	5.93	4.74	7.12	
	μmol/l	8.35	6.68	10.0	Roche Cobas 6000/8000 Online KIMS
	μg/ml	3.99	3.19	4.79	
	μmol/l	10.9	8.72	13.1	Roche Cobas 6000/8000 CEDIA
	μg/ml	5.21	4.17	6.25	
	μmol/l	12.1	9.68	14.5	Abbott Architect (PETINA)
	μg/ml	5.78	4.63	6.93	
	μmol/l	11.7	9.36	14.0	Siemens Atellica CH
	μg/ml	5.59	4.47	6.71	
	μmol/l	12.4	9.92	14.9	Abbott Alinity (CMIA)
	μg/ml	5.93	4.74	7.12	
μmol/l	12.2	9.76	14.6	Abbott Alinity (PETINA)	
μg/ml	5.83	4.67	6.99		

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学	
庆大霉素	μmol/l	11.7	9.36	14.0	Gravimetric	
	μg/ml	5.59	4.47	6.69		
生长激素 (GH)	μU/ml	49.5	39.6	59.4	Siemens Immulite 2000/2500	
	ng/ml	16.5	13.2	19.8		
	μU/ml	51.4	41.1	61.7	Siemens Immulite 1000	
	ng/ml	17.1	13.7	20.5		
	μU/ml	39.0	31.2	46.8	Beckman Dxl800	
	ng/ml	13.0	10.4	15.6		
	μU/ml	42.4	33.9	50.9	Beckman Access	
	ng/ml	14.1	11.3	16.9		
	μU/ml	50.1	40.1	60.1	Tosoh Series	
	ng/ml	16.7	13.4	20.0		
	μU/ml	44.1	35.3	52.9	Diasorin Liaison	
	ng/ml	14.7	11.8	17.6		
	μU/ml	43.4	34.7	52.1	Roche Cobas 4000/E411	
	ng/ml	14.5	11.6	17.4		
	μU/ml	44.6	35.7	53.5	Roche Cobas e601/602	
	ng/ml	14.9	11.9	17.9		
免疫球蛋白 E	μU/ml	42.1	33.7	50.5	Diasorin Liaison XL	
	ng/ml	14.0	11.2	16.8		
	μU/ml	45.3	36.2	54.4	Roche Cobas e402/e801	
	ng/ml	15.1	12.1	18.1		
	μU/ml	45.0	36.0	54.0	Tosoh AIA-CL Series	
	ng/ml	15.0	12.0	18.0		
	免疫球蛋白 E	KIU/l = IU/ml	585	468	702	Abbott Architect
		KIU/l = IU/ml	599	479	719	BioMerieux Vidas
		KIU/l = IU/ml	684	547	821	Siemens Centaur XP/XPT/Classic
		KIU/l = IU/ml	634	507	761	Siemens Immulite 2000/2500
		KIU/l = IU/ml	675	540	810	Siemens Immulite 1000
		KIU/l = IU/ml	674	539	809	Roche Elecsys
		KIU/l = IU/ml	576	461	691	Pharmacia Unicap
		KIU/l = IU/ml	574	459	689	Tosoh Series
		KIU/l = IU/ml	663	530	796	Roche Cobas 4000/E411
		KIU/l = IU/ml	663	530	796	Roche Cobas e601/602
KIU/l = IU/ml		685	548	822	Siemens Centaur CP	
KIU/l = IU/ml		685	548	822	Siemens Atellica IM	
KIU/l = IU/ml		648	518	778	Roche Cobas e402/e801	
KIU/l = IU/ml		576	461	691	Abbott Alinity	
胰岛素	mU/l	20.3	15.2	25.4	Siemens Centaur XP/XPT/Classic	
	μU/ml	20.3	15.2	25.4		
	mU/l	8.25	6.19	10.3	Siemens Immulite 1000	
	μU/ml	8.25	6.19	10.3		
	mU/l	8.84	6.63	11.1	Siemens Immulite 2000	
	μU/ml	8.84	6.63	11.1		
	mU/l	28.6	21.5	35.8	Roche Elecsys	
	μU/ml	28.6	21.5	35.7		

免疫质控品-水平3

货号: IA3111 批号: 2198EC 规格: 12 x 5ml 效期: 2024-12

范围					
分析物	单位	靶值	低值	高值	方法学
胰岛素	mU/l	20.8	15.6	26.0	Beckman Access
	μU/ml	20.8	15.6	26.0	
	mU/l	19.7	14.8	24.6	Abbott Architect
	μU/ml	19.7	14.8	24.6	
	mU/l	36.6	27.5	45.8	Tosoh Series
	μU/ml	36.6	27.5	45.7	
	mU/l	28.3	21.2	35.4	Roche Cobas 4000/E411
	μU/ml	28.3	21.2	35.4	
	mU/l	19.9	14.9	24.9	Beckman Dxl 600/800
	μU/ml	19.9	14.9	24.9	
	mU/l	29.7	22.3	37.1	Diasorin Liaison
	μU/ml	29.7	22.3	37.1	
	mU/l	27.9	20.9	34.9	Roche Cobas 6000/8000
	μU/ml	27.9	20.9	34.9	
	mU/l	25.5	19.1	31.9	Fujirebio Lumipulse G Series
	μU/ml	25.5	19.1	31.9	
	mU/l	21.0	15.8	26.3	Siemens Centaur CP
	μU/ml	21.0	15.8	26.2	
	mU/l	31.7	23.8	39.6	Diasorin Liaison XL
	μU/ml	31.7	23.8	39.6	
mU/l	21.4	16.1	26.8	Siemens Atellica IM	
μU/ml	21.4	16.1	26.7		
mU/l	17.9	13.4	22.4	Abbott Alinity i	
μU/ml	17.9	13.4	22.4		
mU/l	33.2	24.9	41.5	Tosoh AIA-CL Series	
μU/ml	33.2	24.9	41.5		
促黄体生成素(LH)	mU/ml	42.1	33.7	50.5	Abbott Architect
	mU/ml	51.5	41.2	61.8	BioMerieux Vidas
	mU/ml	49.4	39.5	59.3	Siemens Centaur XP/XPT/Classic
	mU/ml	70.6	56.5	84.7	Siemens Immulite 2000
	mU/ml	80.1	64.1	96.1	Siemens Immulite 1000
	mU/ml	41.9	33.5	50.3	Beckman Dxl800
	mU/ml	54.8	43.8	65.8	Roche Elecsys
	mU/ml	63.3	50.6	76.0	Diasorin Liaison
	mU/ml	42.9	34.3	51.5	Beckman Access
	mU/ml	49.1	39.3	58.9	Tosoh Series
	mU/ml	59.5	47.6	71.4	Vitros ECi
	mU/ml	56.3	45.0	67.6	Roche Cobas 4000/E411
	mU/ml	55.6	44.5	66.7	Roche Cobas e601/602
	mU/ml	49.3	39.4	59.2	Fujirebio Lumipulse G Series
	mU/ml	48.0	38.4	57.6	Siemens Centaur CP
	mU/ml	65.3	52.2	78.4	Diasorin Liaison XL
	mU/ml	50.0	40.0	60.0	Siemens Atellica IM
	mU/ml	41.4	33.1	49.7	Abbott Alinity i
	mU/ml	55.4	44.3	66.5	Roche Cobas e402/e801

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
Oestradiol	pmol/l	1063	850	1276	Abbott Architect
	pg/ml	289	231	347	
	pmol/l	2828	2262	3394	BioMerieux Vidas
	pg/ml	769	615	923	
	pmol/l	2055	1644	2466	Siemens Centaur XP/XPT/Classic
	pg/ml	559	447	671	
	pmol/l	2045	1636	2454	Siemens Immulite 2000/2500
	pg/ml	556	445	667	
	pmol/l	1917	1534	2300	Siemens Immulite 1000
	pg/ml	521	417	625	
	pmol/l	1146	917	1375	Beckman Dxl800
	pg/ml	312	249	375	
	pmol/l	1194	955	1433	Roche Elecsys (Estradiol II)
	pg/ml	325	260	390	
	pmol/l	1244	995	1493	Beckman Access
	pg/ml	338	271	405	
	pmol/l	1518	1214	1822	Tosoh (E2)
	pg/ml	413	330	496	
	pmol/l	2408	1926	2890	Vitros ECi
	pg/ml	655	524	786	
	pmol/l	927	742	1112	Diasorin Liaison
	pg/ml	252	202	302	
	pmol/l	1228	982	1474	Roche Cobas e601/602 (Estradiol II)
	pg/ml	334	267	401	
	pmol/l	1237	990	1484	Roche Cobas 4000/E411 (Estradiol II)
	pg/ml	336	269	403	
	pmol/l	2087	1670	2504	Siemens Centaur XP/XPT/Classic Enhanced (eE2)
	pg/ml	568	454	682	
	pmol/l	2107	1686	2528	Siemens Centaur XP/XPT/Classic E26 III
	pg/ml	573	459	687	
	pmol/l	1480	1184	1776	Tosoh (iE2-hsE2)
	pg/ml	403	322	484	
	pmol/l	1775	1420	2130	Siemens Centaur CP
	pg/ml	483	386	580	
	pmol/l	2163	1730	2596	Siemens Centaur CP Enhanced (eE2)
	pg/ml	588	471	705	
	pmol/l	1241	993	1489	Roche Cobas e601/602 (Estradiol III)
	pg/ml	338	270	406	
	pmol/l	1229	983	1475	Roche Elecsys (Estradiol III)
	pg/ml	334	267	401	
	pmol/l	1246	997	1495	Roche Cobas 4000/E411 (Estradiol III)
	pg/ml	339	271	407	
	pmol/l	874	699	1049	Diasorin Liaison XL
	pg/ml	238	190	286	
	pmol/l	1243	994	1492	Beckman Access / Access 2 Sensitive B84493
	pg/ml	338	270	406	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学	
雌二醇 (E2)	pmol/l	1117	894	1340	Beckman DxI600/800 Sensitive B84493	
	pg/ml	304	243	365		
	pmol/l	2171	1737	2605	Siemens Atellica IM	
	pg/ml	591	472	710		
	pmol/l	1037	830	1244	Abbott Alinity i	
	pg/ml	282	226	338		
	pmol/l	1234	987	1481	Roche Cobas e801 (Estradiol III)	
	pg/ml	336	268	404		
	pmol/l	1398	1118	1678	Tosoh AIA-CL Series	
	pg/ml	380	304	456		
	扑热息痛	mmol/l	0.486	0.389	0.583	Abbott Architect
		mg/l	73.5	58.9	88.1	
mmol/l		0.484	0.387	0.581	Roche Cobas c501/502/503/701/702	
mg/l		73.2	58.6	87.8		
mmol/l		0.504	0.403	0.605	Beckman AU Instruments	
mg/l		76.3	61.0	91.6		
mmol/l		0.511	0.409	0.613	Abbott Alinity	
mg/l		77.3	61.9	92.7		
mmol/l		0.488	0.390	0.586	Siemens Atellica CH	
mg/l		73.8	59.0	88.6		
mmol/l		0.529	0.420	0.630	Gravimetric	
mg/l		80.0	63.5	95.3		
苯巴比妥	μmol/l	124	99.0	149	Abbott Architect (CMIA)	
	μg/ml	28.8	23.0	34.6		
	μmol/l	114	91.2	137	Roche Cobas e601/602/801	
	μg/ml	26.4	21.2	31.6		
	μmol/l	123	98.4	148	Siemens Immulite 2000/2500	
	μg/ml	28.5	22.8	34.2		
	μmol/l	121	96.8	145	Beckman AU400/600/640/2700/5400	
	μg/ml	28.1	22.5	33.7		
	μmol/l	173	138	208	Siemens Centaur XP/XPT/Classic	
	μg/ml	40.1	32.0	48.2		
	μmol/l	117	93.6	140	Roche Integra	
	μg/ml	27.1	21.7	32.5		
	μmol/l	118	94.4	142	Roche Cobas C303/501/502/503/701	
	μg/ml	27.4	21.9	32.9		
	μmol/l	114	91.2	137	Abbott Architect (PETINA)	
	μg/ml	26.4	21.2	31.6		
	μmol/l	124	99.0	149	Siemens Atellica CH	
	μg/ml	28.8	23.0	34.6		
μmol/l	108	86.4	130	Abbott Alinity (PETINA)		
μg/ml	25.1	20.0	30.2			
μmol/l	121	96.8	145	Gravimetric		
μg/ml	28.1	22.5	33.6			
苯妥英钠	μmol/l	94.9	75.9	114	Abbott Architect	
	μg/ml	24.0	19.2	28.8		

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
苯妥英钠	μmol/l	119	95.2	143	Siemens Centaur XP/XPT/Classic
	μg/ml	30.0	24.0	36.0	
	μmol/l	105	84.0	126	Siemens Dimension
	μg/ml	26.5	21.2	31.8	
	μmol/l	97.4	77.9	117	Siemens Immulite 2000/2500
	μg/ml	24.6	19.7	29.5	
	μmol/l	91.7	73.4	110	Roche Integra
	μg/ml	23.2	18.5	27.9	
	μmol/l	98.1	78.5	118	Vitros
	μg/ml	24.8	19.8	29.8	
	μmol/l	92.5	74.0	111	Roche Cobas 6000/8000
	μg/ml	23.4	18.7	28.1	
	μmol/l	94.0	75.2	113	Beckman DxC
	μg/ml	23.7	19.0	28.4	
	μmol/l	95.0	76.0	114	Beckman AU Instruments
	μg/ml	24.0	19.2	28.8	
	μmol/l	90.7	72.6	109	EMIT
	μg/ml	22.9	18.3	27.5	
	μmol/l	101	80.8	121	Siemens Atellica CH
	μg/ml	25.5	20.4	30.6	
μmol/l	102	81.6	122	Abbott Alinity	
μg/ml	25.8	20.6	31.0		
μmol/l	99.8	79.8	120	Gravimetric	
μg/ml	25.2	20.1	30.3		
扑米酮	μmol/l	32.7	26.2	39.2	Gravimetric
	μg/ml	7.14	5.72	8.56	
孕酮	nmol/l	87.6	70.1	105	Abbott Architect
	ng/ml	27.4	21.9	32.9	
	nmol/l	134	107	161	BioMerieux Vidas
	ng/ml	41.9	33.5	50.3	
	nmol/l	67.8	54.2	81.4	Siemens Centaur XP/XPT/Classic
	ng/ml	21.2	17.0	25.4	
	nmol/l	71.1	56.9	85.3	Siemens Immulite 2000/2500
	ng/ml	22.3	17.8	26.8	
	nmol/l	80.5	64.4	96.6	Siemens Immulite 1000
	ng/ml	25.2	20.2	30.2	
	nmol/l	86.2	69.0	103	Beckman Dxl Prog Cal 33555
	ng/ml	27.0	21.6	32.4	
	nmol/l	89.2	71.4	107	Roche Elecsys (Prog. II)
	ng/ml	27.9	22.3	33.5	
	nmol/l	87.0	69.6	104	Beckman Access Prog Cal 33555
	ng/ml	27.2	21.8	32.6	
nmol/l	101	80.8	121	Tosoh Series	
ng/ml	31.6	25.3	37.9		
nmol/l	98.9	79.1	119	Vitros ECi	
ng/ml	31.0	24.8	37.2		

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
孕酮	nmol/l	91.5	73.2	110	Roche Cobas 4000/e411 (Prog. II)
	ng/ml	28.6	22.9	34.3	
	nmol/l	92.4	73.9	111	Roche Cobas e601/602 (Prog. II)
	ng/ml	28.9	23.1	34.7	
	nmol/l	80.4	64.3	96.5	Beckman Dxl Prog DE Cal A80773
	ng/ml	25.2	20.1	30.3	
	nmol/l	85.2	68.2	102	Beckman Access Prog DE Cal A80773
	ng/ml	26.7	21.3	32.1	
	nmol/l	64.6	51.7	77.5	Siemens Centaur CP
	ng/ml	20.2	16.2	24.2	
	nmol/l	91.6	73.3	110	Roche Cobas e601/602 (Prog. III)
	ng/ml	28.7	22.9	34.5	
	nmol/l	92.4	73.9	111	Roche Elecsys (Prog. III)
	ng/ml	28.9	23.1	34.7	
	nmol/l	91.4	73.1	110	Roche Cobas 4000/e411 (Prog. III)
	ng/ml	28.6	22.9	34.3	
	nmol/l	81.3	65.0	97.6	Diasorin Liaison XL
	ng/ml	25.4	20.3	30.5	
nmol/l	82.4	65.9	98.9	Abbott Alinity i	
ng/ml	25.8	20.6	31.0		
nmol/l	92.7	74.2	111	Roche Cobas e801 (Prog. III)	
ng/ml	29.0	23.2	34.8		
nmol/l	99.0	79.2	119	Tosoh AIA-CL Series	
ng/ml	31.0	24.8	37.2		
催乳素 (PRL)	μU/ml	1193	954	1432	Abbott Architect
	μU/ml	1363	1090	1636	BioMerieux Vidas
	μU/ml	976	781	1171	Siemens Centaur XP/XPT/Classic
	μU/ml	1068	854	1282	Siemens Immulite 2000/2500
	μU/ml	1180	944	1416	Siemens Immulite 1000
	μU/ml	943	754	1132	Beckman Dxl800
	μU/ml	1088	870	1306	Roche Elecsys
	μU/ml	1072	858	1286	Diasorin Liaison
	μU/ml	978	782	1174	Beckman Access
	μU/ml	1124	899	1349	Tosoh Series
	μU/ml	983	786	1180	Vitros ECi
	μU/ml	1059	847	1271	Roche Elecsys/E170/Cobas 6000/8000 PRLII
	μU/ml	1147	918	1376	Immunotech IRMA
	μU/ml	1107	886	1328	Roche Cobas 4000/E411
	μU/ml	1062	850	1274	Roche Cobas e601/602
	μU/ml	944	755	1133	Siemens Centaur CP
	μU/ml	1178	942	1414	Diasorin Liaison XL
	μU/ml	983	786	1180	Siemens Atellica IM
	μU/ml	1213	970	1456	Abbott Alinity i
	μU/ml	1055	844	1266	Roche Cobas e801 PRL II
游离PSA	ng/ml = μg/l	26.4	19.8	33.0	Abbott Architect
	ng/ml = μg/l	26.7	20.0	33.4	BioMerieux Vidas

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
游离PSA	ng/ml = µg/l	21.0	15.8	26.3	Siemens Centaur XP/XPT/Classic
	ng/ml = µg/l	20.2	15.2	25.3	Siemens Immulite 2000/2500
	ng/ml = µg/l	22.2	16.7	27.8	Siemens Immulite 1000
	ng/ml = µg/l	22.4	16.8	28.0	Beckman DXI standardised to Hybritech
	ng/ml = µg/l	16.6	12.5	20.8	Roche Elecsys
	ng/ml = µg/l	19.7	14.8	24.6	Diasorin Liaison
	ng/ml = µg/l	21.5	16.1	26.9	Beckman Access standardised to Hybritech
	ng/ml = µg/l	16.2	12.2	20.3	Tosoh Series
	ng/ml = µg/l	17.0	12.8	21.3	Roche Cobas 4000/E411
	ng/ml = µg/l	16.9	12.7	21.1	Roche Cobas e601/602
	ng/ml = µg/l	18.0	13.5	22.5	Siemens Dimension Vista
	ng/ml = µg/l	20.9	15.7	26.1	Beckman Dxl800 standardised to WHO IRP96/670
	ng/ml = µg/l	15.5	11.6	19.4	Ortho Vitros 3600/5600/ECi PSA II
	ng/ml = µg/l	18.2	13.7	22.8	Fujirebio Lumipulse G Series
	ng/ml = µg/l	17.6	13.2	22.0	Diasorin Liaison XL
	ng/ml = µg/l	21.3	16.0	26.6	Siemens Atellica IM
ng/ml = µg/l	24.3	18.2	30.4	Abbott Alinity i	
ng/ml = µg/l	16.4	12.3	20.5	Roche Cobas e402/e801	
PSA 总	ng/ml = µg/l	17.7	13.3	22.1	Abbott Architect
	ng/ml = µg/l	25.1	18.8	31.4	BioMerieux Vidas
	ng/ml = µg/l	24.9	18.7	31.1	bioMerieux VIDAS TPSA
	ng/ml = µg/l	22.2	16.7	27.8	Siemens Centaur XP/XPT/Classic
	ng/ml = µg/l	23.6	17.7	29.5	Siemens Immulite 2000/2500 1st generation
	ng/ml = µg/l	23.5	17.6	29.4	Siemens Immulite 2000/2500 3rd generation
	ng/ml = µg/l	22.3	16.7	27.9	Siemens Immulite 1000 1st generation
	ng/ml = µg/l	22.4	16.8	28.0	Beckman DXI standardised to Hybritech
	ng/ml = µg/l	17.2	12.9	21.5	Beckman Dxl800 standardised to WHO IRP96/670
	ng/ml = µg/l	22.6	17.0	28.3	Roche Elecsys
	ng/ml = µg/l	21.8	16.4	27.3	Diasorin Liaison
	ng/ml = µg/l	23.1	17.3	28.9	Beckman Access standardised to Hybritech
	ng/ml = µg/l	22.8	17.1	28.5	Beckman Access standardised to WHO IRP96/670
	ng/ml = µg/l	17.6	13.2	22.0	Tosoh Series
	ng/ml = µg/l	22.9	17.2	28.6	Vitros ECi
	ng/ml = µg/l	24.5	18.4	30.6	CIS RIA Coated tube
	ng/ml = µg/l	23.0	17.3	28.8	Roche Cobas 4000/E411
	ng/ml = µg/l	22.5	16.9	28.1	Roche Cobas e601/602
	ng/ml = µg/l	17.0	12.8	21.3	Fujirebio Lumipulse G Series
	ng/ml = µg/l	22.6	17.0	28.3	Ortho Vitros 3600/5600/ECi PSA II
	ng/ml = µg/l	20.8	15.6	26.0	Siemens Centaur CP
	ng/ml = µg/l	21.6	16.2	27.0	Diasorin Liaison XL
	ng/ml = µg/l	23.7	17.8	29.6	Siemens Atellica IM
	ng/ml = µg/l	17.5	13.1	21.9	Abbott Alinity i
ng/ml = µg/l	22.1	16.6	27.6	Roche Cobas e402/e801	
水杨酸	mmol/l	2.83	2.26	3.40	Abbott Architect
	mg/dl	39.1	31.2	47.0	
	mmol/l	2.82	2.26	3.38	Roche Cobas c501/502/503
	mg/dl	38.9	31.2	46.6	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
水杨酸	mmol/l	2.90	2.32	3.48	EMIT
	mg/dl	40.1	32.0	48.2	
	mmol/l	2.88	2.30	3.46	Gravimetric
	mg/dl	39.8	31.8	47.8	
性激素结合球蛋白 (SHBG)	nmol/l	39.2	31.4	47.0	Siemens Immulite 2000/2500
	nmol/l	41.3	33.0	49.6	Abbott Architect
	nmol/l	42.2	33.8	50.6	Beckman Access
	nmol/l	40.9	32.7	49.1	Beckman Dxl800
	nmol/l	40.7	32.6	48.8	Siemens Atellica IM
	nmol/l	39.4	31.5	47.3	Roche analysers - 413417/15
	nmol/l	37.5	30.0	45.0	Roche Cobas e801 - 413429
	nmol/l	41.2	33.0	49.4	Abbott Alinity i
甲状腺素结合力	%U	38.9	31.9	45.9	Siemens Immulite 2000
睾酮	nmol/l	29.3	23.4	35.2	Abbott Architect
	ng/ml	8.44	6.74	10.1	
	ng/dl	844	674	1014	
	nmol/l	26.8	21.4	32.2	BioMerieux Vidas
	ng/ml	7.72	6.16	9.28	
	ng/dl	772	616	928	
	nmol/l	24.8	19.8	29.8	Siemens Centaur XP/XPT/Classic
	ng/ml	7.14	5.70	8.58	
	ng/dl	714	570	858	
	nmol/l	28.9	23.1	34.7	CIS RIA Coated tube
	ng/ml	8.32	6.65	9.99	
	ng/dl	832	665	999	
	nmol/l	32.1	25.7	38.5	Siemens Immulite 2000/2500
	ng/ml	9.24	7.40	11.1	
	ng/dl	924	740	1108	
	nmol/l	23.4	18.7	28.1	Beckman Dxl800
	ng/ml	6.74	5.39	8.09	
	ng/dl	674	539	809	
	nmol/l	25.9	20.7	31.1	Beckman Access
	ng/ml	7.46	5.96	8.96	
ng/dl	746	596	896		
nmol/l	37.8	30.2	45.4	Tosoh Series	
ng/ml	10.9	8.70	13.1		
ng/dl	1089	870	1308		
nmol/l	38.7	31.0	46.4	Vitros ECi	
ng/ml	11.1	8.93	13.3		
ng/dl	1115	893	1337		
nmol/l	28.4	22.7	34.1	Roche Cobas e601/602 Test II	
ng/ml	8.18	6.54	9.82		
ng/dl	818	654	982		
nmol/l	28.5	22.8	34.2	Roche Elecsys Test II	
ng/ml	8.21	6.57	9.85		
ng/dl	821	657	985		

RANDEX

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
睾酮	nmol/l	28.3	22.6	34.0	Roche Cobas 4000/e411 Test II
	ng/ml	8.15	6.51	9.79	
	ng/dl	815	651	979	
	nmol/l	28.4	22.7	34.1	Abbott Architect 2nd Generation
	ng/ml	8.18	6.54	9.82	
	ng/dl	818	654	982	
	nmol/l	25.8	20.6	31.0	Siemens Centaur CP
	ng/ml	7.43	5.93	8.93	
	ng/dl	743	593	893	
	nmol/l	26.4	21.1	31.7	BioMerieux Vidas Testosterone II
	ng/ml	7.60	6.08	9.12	
	ng/dl	760	608	912	
	nmol/l	24.8	19.8	29.8	Siemens Centaur TST II
	ng/ml	7.14	5.70	8.58	
	ng/dl	714	570	858	
	nmol/l	25.3	20.2	30.4	Siemens Atellica IM
	ng/ml	7.29	5.82	8.76	
	ng/dl	729	582	876	
	nmol/l	28.3	22.6	34.0	Roche Cobas e801 Test II
	ng/ml	8.15	6.51	9.79	
	ng/dl	815	651	979	
nmol/l	34.0	27.2	40.8	Tosoh AIA-CL Series	
ng/ml	9.79	7.83	11.8		
ng/dl	979	783	1175		
nmol/l	29.3	23.4	35.2	Abbott Alinity i	
ng/ml	8.44	6.74	10.1		
ng/dl	844	674	1014		
茶碱	μmol/l	152	122	182	Siemens Centaur XP/XPT/Classic
	μg/ml	27.4	22.0	32.8	
	μmol/l	155	124	186	Siemens Dimension
	μg/ml	27.9	22.3	33.5	
	μmol/l	168	134	202	Roche Integra
	μg/ml	30.3	24.1	36.5	
	μmol/l	146	117	175	Abbott Architect (CMIA)
	μg/ml	26.3	21.1	31.5	
	μmol/l	162	130	194	Roche Cobas 6000/8000
	μg/ml	29.2	23.4	35.0	
	μmol/l	167	134	200	Beckman DxC
	μg/ml	30.1	24.1	36.1	
	μmol/l	142	114	170	Beckman AU400/600/640/2700/5400
	μg/ml	25.6	20.5	30.7	
	μmol/l	159	127	191	Abbott Alinity Enzyme-Immu
μg/ml	28.7	22.9	34.5		
μmol/l	166	133	199	Gravimetric	
μg/ml	29.9	24.0	35.9		

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
甲状腺球蛋白 (TG)	ng/ml	37.3	28.0	46.6	Siemens Immulite 1000
	ng/ml	37.9	28.4	47.4	Siemens Immulite 2000
	ng/ml	72.3	54.2	90.4	Roche Elecsys
	ng/ml	58.4	43.8	73.0	Beckman Dxl800
	ng/ml	70.3	52.7	87.9	Roche Cobas 4000/E411
	ng/ml	68.5	51.4	85.6	Roche Cobas e601/602
	ng/ml	60.5	45.4	75.6	Beckman Access
	ng/ml	80.1	60.1	100	Abbott Architect 2000i
	ng/ml	67.2	50.4	84.0	Tosoh Series
	ng/ml	68.8	51.6	86.0	Roche Cobas e601/602 Tg II
	ng/ml	69.0	51.8	86.3	Roche Elecsys TGII
	ng/ml	69.4	52.1	86.8	Roche Cobas 4000/E411 TGII
	ng/ml	59.4	44.6	74.3	Diasorin Liaison XL
	ng/ml	66.1	49.6	82.6	Roche Cobas e801 Tg II
	ng/ml	60.7	45.5	75.9	DiaSorin Liaison XL TG II Gen
ng/ml	78.0	58.5	97.5	Abbott Alinity i	
促甲状腺激素 (TSH)	μU/ml = mIU/l	18.4	14.7	22.1	Abbott Architect
	μU/ml = mIU/l	19.2	15.4	23.0	bioMerieux VIDAS TSH
	μU/ml = mIU/l	16.9	13.5	20.3	bioMerieux VIDAS TSH3 Ultrasensitive
	μU/ml = mIU/l	23.5	18.8	28.2	Siemens Centaur XP/XPT/Classic
	μU/ml = mIU/l	25.1	20.1	30.1	CIS IRMA
	μU/ml = mIU/l	23.6	18.9	28.3	Siemens Centaur XP/XPT/Classic 3rd Generation
	μU/ml = mIU/l	18.8	15.0	22.6	Siemens Dimension
	μU/ml = mIU/l	21.6	17.3	25.9	Siemens Immulite 2000/2500
	μU/ml = mIU/l	21.4	17.1	25.7	Siemens Immulite 1000
	μU/ml = mIU/l	20.3	16.2	24.4	Beckman Dxl800 Fast TSH
	μU/ml = mIU/l	20.6	16.5	24.7	Beckman Dxl800
	μU/ml = mIU/l	22.1	17.7	26.5	Roche Elecsys
	μU/ml = mIU/l	30.1	24.1	36.1	Diasorin Liaison
	μU/ml = mIU/l	20.0	16.0	24.0	Beckman Access Fast TSH
	μU/ml = mIU/l	20.7	16.6	24.8	Beckman Access hyperTSH 3rd Generation
	μU/ml = mIU/l	24.1	19.3	28.9	Tosoh Series
	μU/ml = mIU/l	28.2	22.6	33.8	Vitros ECI
	μU/ml = mIU/l	22.0	17.6	26.4	Roche Cobas 4000/E411
	μU/ml = mIU/l	22.0	17.6	26.4	Roche Cobas e601/602
	μU/ml = mIU/l	20.8	16.6	25.0	Beckman Dxl800 Hyper TSH
	μU/ml = mIU/l	20.0	16.0	24.0	Siemens Dimension Vista LOCI
	μU/ml = mIU/l	23.9	19.1	28.7	Siemens Centaur XP/XPT/Classic TSH3-Ultra
	μU/ml = mIU/l	22.7	18.2	27.2	Siemens Centaur CP
	μU/ml = mIU/l	22.2	17.8	26.6	Siemens Centaur CP TSH3-Ultra
	μU/ml = mIU/l	20.3	16.2	24.4	Beckman Dxl 600/800 Access (3rd IS)
	μU/ml = mIU/l	32.8	26.2	39.4	Diasorin Liaison XL
	μU/ml = mIU/l	23.6	18.9	28.3	Siemens Atellica IM
μU/ml = mIU/l	18.2	14.6	21.8	Abbott Alinity i	
μU/ml = mIU/l	20.8	16.6	25.0	Roche Cobas e402/e801	
妥布霉素	μmol/l	10.3	8.24	12.4	Gravimetric
	μg/ml	4.82	3.86	5.80	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
总 β-hCG	mU/ml=IU/l	260	208	312	Abbott Architect
	IU/ml	0.260	0.208	0.312	
	mU/ml=IU/l	471	377	565	BioMerieux Vidas
	IU/ml	0.471	0.377	0.565	
	mU/ml=IU/l	227	182	272	Siemens Centaur XP/XPT/Classic
	IU/ml	0.227	0.182	0.272	
	mU/ml=IU/l	456	365	547	Siemens Dimension
	IU/ml	0.456	0.365	0.547	
	mU/ml=IU/l	523	418	628	Siemens Immulite 2000
	IU/ml	0.523	0.418	0.628	
	mU/ml=IU/l	547	438	656	Siemens Immulite 1000
	IU/ml	0.547	0.438	0.656	
	mU/ml=IU/l	268	214	322	Beckman Dxl800
	IU/ml	0.268	0.214	0.322	
	mU/ml=IU/l	243	194	292	Diasorin Liaison
	IU/ml	0.243	0.194	0.292	
	mU/ml=IU/l	305	244	366	Roche Cobas Core EIA
	IU/ml	0.305	0.244	0.366	
	mU/ml=IU/l	260	208	312	Beckman Access
	IU/ml	0.260	0.208	0.312	
	mU/ml=IU/l	471	377	565	Tosoh Series
	IU/ml	0.471	0.377	0.565	
	mU/ml=IU/l	461	369	553	Vitros ECi
	IU/ml	0.461	0.369	0.553	
	mU/ml=IU/l	242	194	290	Roche hCG STAT (Intact)
	IU/ml	0.242	0.194	0.290	
	mU/ml=IU/l	310	248	372	Roche hCG+Beta
	IU/ml	0.310	0.248	0.372	
	mU/ml=IU/l	258	206	310	Beckman Access Total BhCG (5th IS)
	IU/ml	0.258	0.206	0.310	
mU/ml=IU/l	270	216	324	Beckman Dxl Total BhCG (5th IS)	
IU/ml	0.270	0.216	0.324		
mU/ml=IU/l	223	178	268	Siemens Centaur CP	
IU/ml	0.223	0.178	0.268		
mU/ml=IU/l	260	208	312	Diasorin Liaison XL	
IU/ml	0.260	0.208	0.312		
mU/ml=IU/l	243	194	292	Siemens Atellica IM	
IU/ml	0.243	0.194	0.292		
mU/ml=IU/l	260	208	312	Abbott Alinity i	
IU/ml	0.260	0.208	0.312		
mU/ml=IU/l	295	236	354	Roche hCG+Beta e801	
IU/ml	0.295	0.236	0.354		
mU/ml=IU/l	559	447	671	Tosoh AIA-CL Series	
IU/ml	0.559	0.447	0.671		
总 T3	nmol/l	4.69	3.52	5.86	Abbott Architect
	ng/ml	3.05	2.29	3.81	
	ng/dl	305	229	381	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
总 T3	nmol/l	5.22	3.92	6.53	BioMerieux Vidas
	ng/ml	3.40	2.55	4.25	
	ng/dl	340	255	425	
	nmol/l	6.01	4.51	7.51	Siemens Centaur XP/XPT/Classic
	ng/ml	3.91	2.94	4.88	
	ng/dl	391	294	488	
	nmol/l	5.46	4.10	6.83	Siemens Immulite 2000/2500
	ng/ml	3.55	2.67	4.43	
	ng/dl	355	267	443	
	nmol/l	5.31	3.98	6.64	Siemens Immulite 1000
	ng/ml	3.46	2.59	4.33	
	ng/dl	346	259	433	
	nmol/l	5.30	3.98	6.63	Beckman Dxl800
	ng/ml	3.45	2.59	4.31	
	ng/dl	345	259	431	
	nmol/l	5.64	4.23	7.05	Roche Elecsys
	ng/ml	3.67	2.75	4.59	
	ng/dl	367	275	459	
	nmol/l	5.35	4.01	6.69	Beckman Access
	ng/ml	3.48	2.61	4.35	
	ng/dl	348	261	435	
	nmol/l	5.13	3.85	6.41	Tosoh Series
	ng/ml	3.34	2.51	4.17	
	ng/dl	334	251	417	
	nmol/l	6.69	5.02	8.36	Vitros ECI
	ng/ml	4.36	3.27	5.45	
	ng/dl	436	327	545	
	nmol/l	4.72	3.54	5.90	CIS RIA Coated tube
	ng/ml	3.07	2.30	3.84	
	ng/dl	307	230	384	
nmol/l	5.58	4.19	6.98	Roche Cobas 4000/E411	
ng/ml	3.63	2.73	4.53		
ng/dl	363	273	453		
nmol/l	4.13	3.10	5.16	Immunotech RIA	
ng/ml	2.69	2.02	3.36		
ng/dl	269	202	336		
nmol/l	5.47	4.10	6.84	Roche Cobas e601/602	
ng/ml	3.56	2.67	4.45		
ng/dl	356	267	445		
nmol/l	6.31	4.73	7.89	Siemens Centaur CP	
ng/ml	4.11	3.08	5.14		
ng/dl	411	308	514		
nmol/l	4.38	3.29	5.48	Diasorin Liaison XL	
ng/ml	2.85	2.14	3.56		
ng/dl	285	214	356		

RANDOX

免疫质控品-水平3

货号: IA3111 批号: 2198EC 规格: 12 x 5ml 效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学	
总 T3	nmol/l	4.79	3.59	5.99	Abbott Alinity i	
	ng/ml	3.12	2.34	3.90		
	ng/dl	312	234	390		
	nmol/l	5.89	4.42	7.36	Siemens Atellica IM	
		ng/ml	3.83	2.88		4.78
		ng/dl	383	288		478
	nmol/l	5.91	4.43	7.39	Roche Cobas e402/e801	
		ng/ml	3.85	2.88		4.82
		ng/dl	385	288		482
	nmol/l	5.32	3.99	6.65	Tosoh AIA-CL Series	
		ng/ml	3.46	2.60		4.32
		ng/dl	346	260		432
总 T4	nmol/l	188	141	235	Abbott Architect	
	µg/dl	14.7	11.0	18.4		
	ng/ml	147	110	184		
	nmol/l	206	155	258	BioMerieux Vidas	
		µg/dl	16.1	12.1		20.1
		ng/ml	161	121		201
	nmol/l	210	158	263	Siemens Centaur XP/XPT/Classic	
		µg/dl	16.4	12.3		20.5
		ng/ml	164	123		205
	nmol/l	186	140	233	Siemens Immulite 2000/2500	
		µg/dl	14.5	10.9		18.1
		ng/ml	145	109		181
	nmol/l	188	141	235	Siemens Immulite 1000	
		µg/dl	14.7	11.0		18.4
		ng/ml	147	110		184
	nmol/l	206	155	258	Beckman Dxl800	
		µg/dl	16.1	12.1		20.1
		ng/ml	161	121		201
	nmol/l	175	131	219	Roche Elecsys	
		µg/dl	13.7	10.2		17.2
		ng/ml	137	102		172
	nmol/l	200	150	250	Beckman Access	
		µg/dl	15.6	11.7		19.5
		ng/ml	156	117		195
nmol/l	191	143	239	Tosoh Series		
	µg/dl	14.9	11.2		18.6	
	ng/ml	149	112		186	
nmol/l	202	152	253	Vitros ECi		
	µg/dl	15.8	11.9		19.7	
	ng/ml	158	119		197	
nmol/l	236	177	295	CIS RIA Coated tube		
	µg/dl	18.4	13.8		23.0	
	ng/ml	184	138		230	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
总 T4	nmol/l	173	130	216	Roche Cobas 4000/E411
	µg/dl	13.5	10.1	16.9	
	ng/ml	135	101	169	
	nmol/l	171	128	214	Roche Cobas e601/602
	µg/dl	13.3	9.98	16.6	
	ng/ml	133	100	166	
	nmol/l	147	110	184	Immunotech IRMA
	µg/dl	11.5	8.58	14.4	
	ng/ml	115	85.8	144	
	nmol/l	214	161	268	Siemens Centaur CP
	µg/dl	16.7	12.6	20.8	
	ng/ml	167	126	208	
	nmol/l	182	137	228	Diasorin Liaison XL
	µg/dl	14.2	10.7	17.7	
	ng/ml	142	107	177	
	nmol/l	186	140	233	Abbott Alinity i
	µg/dl	14.5	10.9	18.1	
	ng/ml	145	109	181	
	nmol/l	224	168	280	Siemens Atellica IM
	µg/dl	17.5	13.1	21.9	
	ng/ml	175	131	219	
	nmol/l	169	127	211	Roche Cobas e402/e801
	µg/dl	13.2	9.91	16.5	
	ng/ml	132	99.1	165	
nmol/l	215	161	269	Tosoh AIA-CL Series	
µg/dl	16.8	12.6	21.0		
ng/ml	168	126	210		
丙戊酸	µmol/l	940	752	1128	Abbott Architect (CMIA)
	µg/ml	136	109	163	
	µmol/l	992	794	1190	Siemens Centaur XP/XPT/Classic
	µg/ml	143	115	171	
	µmol/l	897	718	1076	Siemens Dimension
	µg/ml	129	104	154	
	µmol/l	998	798	1198	Siemens Immulite 2000/2500
	µg/ml	144	115	173	
	µmol/l	1031	825	1237	Siemens Immulite 1000
	µg/ml	149	119	179	
	µmol/l	934	747	1121	Roche Integra
	µg/ml	135	108	162	
	µmol/l	974	779	1169	Vitros
	µg/ml	141	112	170	
	µmol/l	938	750	1126	Roche Cobas 6000/8000
	µg/ml	135	108	162	
	µmol/l	968	774	1162	Beckman DxC
	µg/ml	140	112	168	
	µmol/l	966	773	1159	Beckman AU400/600/640/2700/5400
	µg/ml	139	112	166	

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学	
丙戊酸	μmol/l	1013	810	1216	Siemens ADVIA 1200/1650/1800/2400	
	μg/ml	146	117	175		
	μmol/l	921	737	1105	Abbott Architect (PETINA)	
	μg/ml	133	106	160		
	μmol/l	996	797	1195	EMIT	
	μg/ml	144	115	173		
	μmol/l	913	730	1096	Siemens Atellica CH	
	μg/ml	132	105	159		
	μmol/l	904	723	1085	Abbott Alinity (CMIA)	
	μg/ml	130	104	156		
	μmol/l	902	722	1082	Abbott Alinity (PETINA)	
	μg/ml	130	104	156		
	μmol/l	1116	893	1339	Gravimetric	
	μg/ml	161	129	193		
	万古霉素	μmol/l	10.7	8.56	12.8	Roche Integra
		μg/ml	15.9	12.7	19.1	
μmol/l		9.64	7.71	11.6	Ortho Vitros Microslide Systems	
μg/ml		14.3	11.5	17.1		
μmol/l		12.4	9.92	14.9	Abbott Architect (CMIA)	
μg/ml		18.4	14.7	22.1		
μmol/l		9.76	7.81	11.7	Siemens Centaur XP/XPT/Classic	
μg/ml		14.5	11.6	17.4		
μmol/l		11.1	8.88	13.3	Siemens Dimension	
μg/ml		16.5	13.2	19.8		
μmol/l		10.9	8.72	13.1	Roche Cobas 6000/8000	
μg/ml		16.2	13.0	19.4		
μmol/l		9.49	7.59	11.4	Beckman DxC	
μg/ml		14.1	11.3	16.9		
μmol/l		10.6	8.48	12.7	Thermo Scientific QMS Vancomycin	
μg/ml		15.7	12.6	18.8		
μmol/l		10.7	8.56	12.8	Abbott Architect (PETINA)	
μg/ml		15.9	12.7	19.1		
μmol/l		11.1	8.88	13.3	Siemens Atellica CH	
μg/ml		16.5	13.2	19.8		
μmol/l	11.0	8.80	13.2	Abbott Alinity (CMIA)		
μg/ml	16.3	13.1	19.5			
μmol/l	10.8	8.64	13.0	Abbott Alinity (PETINA)		
μg/ml	16.0	12.8	19.2			
μmol/l	10.8	8.64	13.0	Gravimetric		
μg/ml	16.0	12.8	19.3			
维生素 B12	pmol/l	515	412	618	Abbott Architect	
	pg/ml	698	558	838		
	pmol/l	625	500	750	Siemens Centaur XP/XPT/Classic	
	pg/ml	847	678	1016		
	pmol/l	589	471	707	Siemens Immulite 2000	
pg/ml	798	638	958			

免疫质控品-水平3

货号: IA3111

批号: 2198EC

规格: 12 x 5ml

效期: 2024-12

范围

分析物	单位	靶值	低值	高值	方法学
维生素 B12	pmol/l	489	391	587	Beckman Dxl800
	pg/ml	663	530	796	
	pmol/l	629	503	755	Roche Elecsys
	pg/ml	852	682	1022	
	pmol/l	514	411	617	Beckman Access
	pg/ml	696	557	835	
	pmol/l	762	610	914	Tosoh Series
	pg/ml	1033	827	1239	
	pmol/l	638	510	766	Vitros ECI
	pg/ml	864	691	1037	
	pmol/l	582	466	698	Siemens Dimension
	pg/ml	789	631	947	
	pmol/l	707	566	848	Siemens Centaur CP
	pg/ml	958	767	1149	
	pmol/l	634	507	761	Roche Cobas 4000/e411 Gen II
	pg/ml	859	687	1031	
	pmol/l	653	522	784	Roche Cobas e601/602 Gen II
	pg/ml	885	707	1063	
	pmol/l	537	430	644	Diasorin Liaison XL
	pg/ml	728	583	873	
pmol/l	665	532	798	Siemens Atellica IM	
pg/ml	901	721	1081		
pmol/l	464	371	557	Abbott Alinity i	
pg/ml	629	503	755		
pmol/l	636	509	763	Roche Cobas e801 Gen II	
pg/ml	862	690	1034		