

LIQUID BNP CONTROL (BECKMAN ACCESS / BECKMAN DxI) - LEVEL I (BNP CONTROL I)

CAT. NO. CQ5133
SIZE: 3 x 1 ml

LOT NO. 4030CK
EXPIRY: 2018-03-28

GTIN: 05055273207477

INTENDED USE

This product is intended for *in vitro* diagnostic use, in the quality control of BNP on the Beckman Access and Beckman DxI analysers.

DEVICE DESCRIPTION

The Liquid BNP Controls are supplied at 3 levels, 1, 2 and 3. Target values and ranges are supplied for BNP only.

SAFETY PRECAUTIONS AND WARNINGS

For *in vitro* diagnostic use only. Do not pipette by mouth. Exercise the normal precautions required for handling laboratory reagents.

Avoid ingestion or contact with skin or mucous membranes. In case of skin contact, flush affected area with copious amounts of water. In case of contact with eyes, or if ingested, seek immediate medical attention.

Human source material, from which this product has been derived, has been tested at donor level for the Human Immunodeficiency Virus (HIV 1, HIV 2) antibody, Hepatitis B Surface Antigen (HbsAg), and Hepatitis C Virus (HCV) antibody and found to be NON-REACTIVE. FDA approved methods have been used to conduct these tests. However, since no method can offer complete assurance as to the absence of infectious agents, this material and all patient samples should be handled as though capable of transmitting infectious diseases and disposed of accordingly.

Health and Safety Data Sheets are available on request.

STORAGE AND STABILITY

UNOPENED: Store at +2°C to +8°C. Stable to expiration date printed on individual vials.

OPENED: Store refrigerated (+2°C to +8°C). Liquid BNP Controls are stable for 30 days at +2°C to +8°C, if kept capped in original container and free from contamination. Only the required amount of product should be removed. After use, any residual product should NOT BE RETURNED to the original vial.

PREPARATION FOR USE

The Liquid BNP Controls are supplied ready to use.

MATERIALS PROVIDED

Liquid BNP Control - Level I 3 x 1 ml

MATERIALS REQUIRED BUT NOT PROVIDED

Not applicable.

ASSIGNED VALUES

Each batch of Liquid BNP Controls is submitted to a number of external laboratories. Values are assigned from consensus of results obtained by these laboratories and internal testing conducted at Randox Laboratories Ltd.

ANALYSER	UNITS	TARGET	RANGE
Beckman Access	pg/ml	86.9	69.5 - 104
Beckman Access	pmol/l	25.2	20.2 - 30.2
Beckman DxI	pg/ml	92.4	73.9 - 111
Beckman DxI	pmol/l	26.8	21.4 - 32.2

16 Nov '16 ne

THIS PAGE IS INTENTIONALLY BLANK