

孕妇多项检测用质控品-水平3

货号 MSS5026 批号 7703MS

规格 3 x 1ml 效期 2023-01-28

分析物	单位	靶值	范围		方法学
			低值	高值	
甲胎蛋白	KIU/l = IU/ml	107	85.6	128	Siemens Immulite 1000
	ng/ml	129	104	154	
	KIU/l = IU/ml	102	81.6	122	Perkin Elmer Delfia/Delfia Xpress/AutoDelfia
	ng/ml	123	98.7	147	
	KIU/l = IU/ml	102	81.6	122	Abbott Architect
	ng/ml	123	98.7	147	
	KIU/l = IU/ml	106	84.8	127	Siemens Immulite 2000/2500
	ng/ml	128	103	153	
	KIU/l = IU/ml	115	92.0	138	Siemens Centaur XP/XPT/Classic
	ng/ml	139	111	167	
	KIU/l = IU/ml	95.2	76.2	114	Beckman Dxl800
	ng/ml	115	92.2	138	
KIU/l = IU/ml	91.5	73.2	110	Brahms Kryptor	
ng/ml	111	88.6	133		
KIU/l = IU/ml	100	80.0	120	Beckman Access	
ng/ml	121	96.8	145		
KIU/l = IU/ml	112	89.6	134	Roche Cobas 6000/8000	
ng/ml	136	108	164		
KIU/l = IU/ml	118	94.4	142	Roche Cobas E411	
ng/ml	143	114	172		
游离 -hCG	mU/mL=U/L	194	146	243	Perkin Elmer Delfia/Delfia Xpress/AutoDelfia
	mU/mL=U/L	235	176	294	Siemens Immulite 2000/2500
	mU/mL=U/L	218	164	273	Siemens Immulite 1000
	mU/mL=U/L	198	149	248	Brahms Kryptor
	mU/mL=U/L	171	128	214	Roche Cobas 6000/8000
	mU/mL=U/L	179	134	224	Roche Cobas 4000/E411
抑制素 A	ng/L = pg/ml	668	501	835	Perkin Elmer Delfia/Delfia Xpress/AutoDelfia
	ng/L = pg/ml	417	313	521	Beckman Dxl 600/800
	ng/L = pg/ml	419	314	524	Beckman Access/LXi725
PAPP-A	U/L=mIU/ml	4.12	3.09	5.15	Perkin Elmer Delfia/Delfia Xpress/AutoDelfia
	ng/ml	18559	13919	23199	
	U/L=mIU/ml	12.7	9.53	15.9	Siemens Immulite 1000
	ng/ml	57207	42928	71486	
	U/L=mIU/ml	13.8	10.4	17.3	Siemens Immulite 2000/2500
	ng/ml	62162	46847	77477	
	U/L=mIU/ml	0.304	0.228	0.380	Beckman Dxl 600/800
	ng/ml	1369	1027	1711	
	U/L=mIU/ml	3.59	2.69	4.49	Brahms Kryptor
	ng/ml	16171	12117	20225	
	U/L=mIU/ml	5.43	4.07	6.79	Roche Cobas 6000/8000
	ng/ml	24459	18333	30585	
U/L=mIU/ml	5.40	4.05	6.75	Roche Cobas E411	
ng/ml	24324	18243	30405		

孕妇多项检测用质控品-水平3

货号 MSS5026 批号 7703MS

规格 3 x 1ml 效期 2023-01-28

分析物	单位	靶值	范围		方法学
			低值	高值	
总 -hCG	mU/ml=IU/l	24225	19380	29070	Abbott Architect I2000/8200
	IU/ml	24.2	19.4	29.0	
	mU/ml=IU/l	21946	17557	26335	Siemens Centaur XP/XPT/Classic
	IU/ml	21.9	17.6	26.2	
	mU/ml=IU/l	48071	38457	57685	Perkin Elmer Delfia/Delfia Xpress/AutoDelfia
	IU/ml	48.1	38.5	57.7	
	mU/ml=IU/l	44251	35401	53101	Siemens Immulite 2000/2500
	IU/ml	44.3	35.4	53.2	
	mU/ml=IU/l	51478	41182	61774	Siemens Immulite 1000
	IU/ml	51.5	41.2	61.8	
	mU/ml=IU/l	29580	23664	35496	Beckman Dxl 600/800
	IU/ml	29.6	23.7	35.5	
	mU/ml=IU/l	28767	23014	34520	Beckman Access/LXi725
	IU/ml	28.8	23.0	34.6	
	mU/ml=IU/l	29431	23545	35317	Roche Cobas 6000/8000
	IU/ml	29.4	23.5	35.3	
mU/ml=IU/l	30042	24034	36050	Beckman Access Total BhCG (5th IS)	
IU/ml	30.0	24.0	36.0		
mU/ml=IU/l	31887	25510	38264	Beckman Dxl Total BhCG (5th IS)	
IU/ml	31.9	25.5	38.3		
未结合雌三醇	nmol/L	40.2	30.2	50.3	Perkin Elmer Delfia/Delfia Xpress/AutoDelfia
	ng/mL	11.6	8.71	14.5	
	nmol/L	17.4	13.1	21.8	Siemens Immulite 2000/2500
	ng/mL	5.02	3.78	6.26	
	nmol/L	18.7	14.0	23.4	Siemens Immulite 1000
	ng/mL	5.39	4.04	6.74	
	nmol/L	14.0	10.5	17.5	Beckman Dxl 600/800
	ng/mL	4.04	3.03	5.05	
	nmol/L	8.64	6.48	10.8	Brahms Kryptor
	ng/mL	2.49	1.87	3.11	
	nmol/L	14.4	10.8	18.0	Beckman Access/LXi725
	ng/mL	4.15	3.11	5.19	