

【英文名称】: **LIPID CONTROL - LEVEL 3 (LPD CONTROL 3)**

【中文名称】: 脂类多项检测用质控品（水平 3）

【货号】 (CAT. NUMBER): **LE2670**

【批号】 (LOT NUMBER): **3102CH**

【包装规格】 (SIZE): **5 x 1 ml**

【到期】 (EXPIRY): **2026-10-28**

本产品的准备和使用详见试剂盒内说明书。

Full instructions for the use and preparation of this product are included with the kit.

质控品靶值和范围如下所示。

Target values and ranges are listed below.

ZH: 19 Jun '24 me

LIPID CONTROL LEVEL 3 (LPD CONTROL 3) 脂类多项检测用质控品 (水平3)							
批号 (Lot. No): 3102CH 货号 (Cat. No): LE2670							
包装规格 (Size): 5 x 1 ml 到期 (Expiry): 2026-10-28							
Analyte	分析项目	unit 单位	Target 靶值	Range		方法	methods
				low 低值	high 高值		
Apolipoprotein A-1	载脂蛋白A-1	g/l	1.96	1.61	2.31	免疫比浊法	Immunoturbidimetric
		mg/dl	196	161	231		
		g/l	2.01	1.65	2.37	散射比浊法	Nephelometric
		mg/dl	201	165	237		
Apolipoprotein B	载脂蛋白 B	g/l	1.51	1.24	1.78	免疫比浊法	Immunoturbidimetric
		mg/dl	151	124	178		
		g/l	1.45	1.19	1.71	散射比浊法	Nephelometric
		mg/dl	145	119	171		
Cholesterol	胆固醇	mmol/l	7.49	6.52	8.46	胆固醇氧化酶 - Abell Kendall	Cholesterol Oxidase - Abell Kendall
		mg/dl	289	252	326		
		mmol/l	7.27	6.32	8.22	Siemens Dimension	Siemens Dimension
		mg/dl	281	244	318		
		mmol/l	7.59	6.60	8.58	胆固醇氧化酶 - IDMS	Cholesterol Oxidase - IDMS
		mg/dl	293	255	331		
HDL - Cholesterol	高密度脂蛋白胆固醇 (HDL-C)	mmol/l	2.06	1.75	2.37	直接清除法	Direct Clearance Method
		mg/dl	79.5	67.6	91.4		
		mmol/l	1.78	1.51	2.05	直接高密度脂蛋白免疫分离法	Direct HDL Immunoseparation
		mg/dl	68.7	58.3	79.1		
		mmol/l	1.71	1.45	1.97	直接高密度脂蛋白PEGME法	Direct HDL PEGME
		mg/dl	66.0	56.0	76.0		
		mmol/l	2.10	1.79	2.42	直接高密度脂蛋白PPD法	Direct HDL PPD
		mg/dl	81.1	69.1	93.1		
LDL - Cholesterol	低密度脂蛋白 (LDL)	mmol/l	1.65	1.40	1.90	直接高密度脂蛋白罗氏第4代	Direct HDL Roche 4th Generation
		mg/dl	63.7	54.0	73.4		
		mmol/l	5.37	4.56	6.18	直接清除法	Direct Clearance Method
		mg/dl	207	176	238		
		mmol/l	3.56	3.03	4.09	计算	Calculated
		mg/dl	137	117	157		
Lipoprotein (a)	脂蛋白 (a)	mmol/l	4.28	3.64	4.92	选择性化学沉淀法	Selective detergent methods
		mg/dl	165	141	189		
		mg/dl	25.7	20.6	30.8	免疫比浊法	Immunoturbidimetric
		nmol/l	49.5	39.6	59.4		
Triglycerides	甘油三酸酯	mmol/l	4.14	3.48	4.80	脂肪酶/GPO-PAP (无校正)	Lipase/GPO-PAP no correction
		mg/dl	366	308	424		
		mmol/l	4.17	3.50	4.84	脂肪酶/GK UV (无校正)	Lipase/GK UV no correction
		mg/dl	369	310	428		
		mmol/l	3.85	3.23	4.47	脂肪酶/GPO-PAP 0.11mmol/l 校正	Lipase/GPO-PAP 0.11mmol/l correction
		mg/dl	341	286	396		
		mmol/l	4.27	3.59	4.95	脂肪酶/甘油脱氢酶	Lipase/Glycerol Dehydrogenase
		mg/dl	378	318	438		
		mmol/l	4.11	3.45	4.77	Siemens Dimension	Siemens Dimension
		mg/dl	364	305	423		